

Uro Quick

Urine Analyzer



Features:

- Q High-luminance and white LED, features in long life and good stability.
- Q Big LCD to display abundant contents, and French and English interface are optional.
- Q User-friendly interface to achieve man-machine conversation.
- Q International unit, conventional unit and Unit Range are optional.
- Q Three modes (One-Step testing, Slowly-Series testing and Quickly-Series testing) are optional, which is convenient for different user.
- Q Monitoring the whole testing process, auto-character and audible-control prompt.
- Q Compatible with 8, 10 and 11 items of testing paper.
- Q Standard RS-232 and USB interface.
- Q Built-in thermal printer.

Uro Quick

Urine Analyzer

Introduction

Uroquick Urine Analyzer is a high-precision, intellectual instrument which is researched and developed basing on modern optics, electronics, computer science and other advanced technologies for clinical inspection of urine. And it integrates the advantages of easy and quick operation, exact result and good repeatability.

Technical Specifications

- Q **Testing items** : GLU, BIL, SG, KET, BLD, PRO, URO, NIT, LEU, VC, PH.
- Q **Testing theory** : RGB tricolor testing theory.
- Q **Repetition** : CV < 1%
- Q **Stability** : CV < 1%
- Q **Display** : 2.8" color LCD.
- Q **Workig mode** : One-Step, Slowly – Series and Quickly-Series testing mode.
- Q **Testing speed** : 120 pieces / hour or 60 pieces / hour.
- Q **Data storage** : Storage of 1000 test items, freely query in term of testing time and sample number.
- Q **Printer** : Built-in high-speed thermal printer.
- Q **Interface** : Standard RS-232 intercommunication interface.
- Q **Power supply** : 100~240V, 50/60Hz.
- Q **Dimension** : 240mm(L) x 220mm (W) x 130mm (H).
- Q **Weight** : about 1.8kg.



G-1, Plot # 30, Cauvery Street, Kubera Nagar, Ramnagar South Ext,
Pallikaranai, Chennai - 600 100. Ph : 044 64545999, +91 9884635363
Email: sales@quicklab.in, service@quicklab.in Web :www.quicklab.in

Our Local Distributor is :