

The 5th generation of
Homogeneous CLIA

QUICKLAB[®]
SERVICES PRIVATE LIMITED

EASY OPERATION



Adding sample
to reagent



Reagent
placement



Result comes in
5 mins

EXTENSIVE SCENARIO



Veterinary
hospital



Animal disease
control



Animal breeding
facility



Independent
clinical laboratory



Pet
beauty salon



Customs and
immigration inspection

Dry Micro Chemiluminescence Immunoassay Analyzer



QUICKLAB[®]
SERVICES
Private Limited
SOURCE & SOLUTION FOR YOUR CLINICAL LAB NEEDS.



Plot No 31B, Kavery Street, Ram Nagar South, Madipakkam,
Chennai - 600091, Tamil Nadu, India



info@quicklab.in, sales1, sales2@quicklab.in, service@quicklab.in



+917449092777, +917449093777, +917449094777, +917449095777,



SPECIFICATION

Test time	5 mins	Mesurement principle	Homogeneous chemiluminescence
Capacity	7 test positions	Pipetting volume	10-100μL
Sample	Serum/Plasm/Whole blood	Storage	Up to 10000 results
Precision	CV<3%	Dimensions	325mm x 231mm x 213mm
Interface	7" LCD touchscreen	Weight	≈ 6kg



Lyophilized Beads

- Low transportation cost
- Long shelf-life

Quality assurance

No cold chain


Individual packaging


High sensitive
High precision


Extensiveness on
sample and assay


No pre-treatment
required

ASSAY LIST

ITEM	SAMPLE	ITEM	SAMPLE
c/f Cor	Serum/Plasma	cNT-proBNP	Serum/Plasma
c/f Prog	Serum/Plasma	fNT-proBNP	Serum/Plasma
c/f HbA1c	Whole blood	cD-Dimer	Plasma
c/f 25-OH VD	Serum/Plasma Whole blood	fD-Dimer	Plasma
c/f FT4	Serum/Plasma	cCRP	Serum/Plasma Whole blood
c/f TT4	Serum/Plasma	fSAA	Serum/Plasma
cTSH	Serum/Plasma	cPL	Serum/Plasma
fST2	Serum/Plasma	fPL	Serum/Plasma

c for Canine, f for Feline